

SPECIFIC TEST PARAMETERS - GENERAL

Specific Test Parameters										
General	LIH	ISE	Range							
Test Name:	SAL	<	>	Type:	Others	Operation:	Yes			
Sample:	Volume	50	µL	Dilution	0	µL	Pre-Dilution Rate:	1		
Reagents:	R1 Volume	140	µL	Dilution	0	µL	Min OD		Max OD	
	R2 Volume	0	µL	Dilution	0	µL	L	-2,0	H	2,5
Wavelength:	Pri.	450		Sec.	520		Reagent OD limit:			
							First L	-2,0	First H	2,5
Method:		END					Last L	-2,0	Last H	2,5
Reaktion slope:		+					Dynamic Range:			
Measuring Point 1:	First	0		Last	27		L	11	H	85
Measuring Point 2:	First			Last			Correlation Factor:			
Linearity:			%				A	1	B	0
No Lag Time:							On-board stability period:	90		

Note:

= User defined

The Test Name "SAL " is suggested as an abbreviation for saliva quantification

Reagent: Aqua dest.

SPECIFIC TEST PARAMETERS - RANGE

Specific Test Parameters									
General	LIH	ISE	Range						
Test Name:	SAL	<	>	Type:	Others				
Value/Flag:	Value	Level L:	#	Level H:	#				
Normal Ranges:	Sex	Year	Month	Year	Month	L	H		
□ 1.	#	#	#	#	#	#	#		
□ 2.	#	#	#	#	#	#	#		
□ 3.	#	#	#	#	#	#	#		
□ 4.	#	#	#	#	#	#	#		
□ 5.	#	#	#	#	#	#	#		
□ 6.	#	#	#	#	#	#	#		
7.	None Selected					11	85		
8.	Out of Range					L	H	-9999	9999
Panic Value:	10		86		Unit:	U/ml	Decimal places:	1	

Note:

= User defined

CALIBRATION SPECIFIC - GENERAL

Calibration Specific					
General	ISE				
Test Name:	SAL	<	>	Type:	Others
Calibration Type:	5AB	Formula:	POLYGONAL	Counts:	2
		Process:	Conc		
	Cal. No.	OD	CONC	Factor/OD-L	Factor/OD-H
Point 1:	#		11	-2	2,5
Point 2:	#		27	-2	2,5
Point 3:	#		46	-2	2,5
Point 4:	#		66	-2	2,5
Point 5:	#		85	-2	2,5
Point 6:					
Point 7:					
1-Point Cal. Point:		<input type="checkbox"/>	With CONC-0	Slope Check:	None
MB Type Factor:			Calibration Stability Period:	30	
				Advanced Cal	No

Note:

= User defined

Calibrators: Greiner Bio-One Saliva Quantification Kit, Saliva Calibrators 1 - 5

QUALITY CONTROL SPECIFIC - PRESETTING

QC Specific						
Presetting	Cumulative					
Test Name	SAL		<	>	Type	Others
Control No	Multi / Single	Mean	Standard deviation	Range	Name of Control	
1.	# Single	30,00	0,70	4	Sal Quantification Kit Control 1	
2.	# Single	70,00	0,70	4	Sal Quantification Kit Control 2	
3.						
4.						
5.						
6.						

Note:

= User defined

Controls: Greiner Bio-One Saliva Quantification Kit, Saliva Control 1 and Saliva Control 2